

Laboratory Technical Bulletin

Date:	June 29, 2016	Effective	Date:	June 29, 2016
Title: I	Metamizole (Dipyrone) interference with Cholest	erol, Trigl	yceride,	Uric Acid, Lactic Acid, and
Explanation : Metamizole has been identified as an interferent which may cause false low results in the following assays: Cholesterol, Triglyceride, Uric Acid, Lactic Acid, and Lipase. Venipuncture should be performed before administration of this drug.				
Dr. Dan	iel Sudilovsky, Laboratory Medical Director	Date	61	29/16
Toni Bu	rger, Laboratory Administrative Director	Date	6-2	9-16