



Laboratory Technical Bulletin

Date: April 4, 2018


Effective Date: April 4, 2018

Title: Effective immediately, Osmolality testing is being performed in house at CMC

Test: Serum Osmolality/Urine Osmolality

	Old	New
<i>Perform Site</i>	Mayo	CMC
<i>Method</i>	Freezing Point Depression	Freezing Point Depression
<i>Test Code</i>	OSM, RUOSM	OSM, RUOSM
<i>Reference Range</i>	Serum (275 – 295) mOsm/kg Urine (150 – 1150) mOsm/kg	Serum (275 – 295) mOsm/kg Urine (150 – 1150) mOsm/kg
<i>Specimen</i>	Serum, urine	Serum, urine
<i>Volume Required</i>	1 mL	1 mL
<i>Turnaround time</i>	24-48 hrs	1 hr


 Date 4/4/18
 Dr. Daniel Sudilovsky, Laboratory Medical Director


 Date 4-4-18
 Toni Burger, Laboratory Administrative Director