


Hemosure Collection Guide

**CLEARLY LABEL SPECIMEN WITH PATIENT'S FULL NAME,
DATE OF BIRTH, DATE AND TIME OF COLLECTION.
*Specimen may be rejected if this crucial step is omitted!***

1. Obtain unopened Hemosure collection tube, specimen bag and stool collection container (if available) from healthcare provider.
2. Use any clean dry container or plastic wrap to collect stool sample.
3. Unscrew the purple cap from the sample collection tube.
DO NOT POUR OUT THE LIQUID.
4. Poke spiral applicator into stool at 6 different sites. Use only enough fecal material to cover the tip of the applicator. DO NOT CLUMP, SCOOP, OR FILL THE TUBE.
5. Screw the applicator back into the tube and secure tightly.
6. Fill out personal information on attached label on tubes.
7. Place specimen in specimen bag for transport.
8. Wash hands thoroughly with soap and water.
9. Deliver hemosure collection tube to healthcare professional.
10. Hemosure tubes containing fecal material may be stored 6 days room temperature or 30 days refrigerated.



Attention!

Fill out personal information
on attached label on tubes.

