



Date: 5/10/2021


Effective Date: 5/14/2021


Title: Fetal Maternal Hemorrhage Rapid Screen, Immucor

Test: Fetal Screen

Explanation of change: The fetal screen test will now be performed on different kit manufactured by Immucor. This Kit requires less incubation time and is incubated at room temperature instead of at 37°C. It is more sensitive and is easier for laboratory technologists to interpret. Providers and nursing units will see no difference in reporting.

	<b>Old</b>	<b>New</b>
Perform Site	CMC Blood Bank	CMC Blood Bank
Method	Manual	Manual
Test Code	FS	FS
Reference Range	Qualitative	No change
Specimen	Pink EDTA	Pink EDTA
Volume Required	6.0mL	6.0mL
Turnaround time	Day shift test only	No change


  
 Dr. Elizabeth Plocharczyk, Laboratory Medical Director
   
 Date: 05/11/21


  
 Patricia Lindquist, Assistant Technical Director
   
 Date: 6/11/21