

**Date:** 8/20/2024

**Effective Date:** 8/20/2024


**Title:** INR Reference Range Change

**Test:** PT/INR

**Method:** ACL TOP 550

**Explanation of change:** Effective 8/20/2024, due to a new lot of PT reagent and a normal patient study, the INR reference range has changed. The PT reference range will remain the same.

	<b>Current</b>	<b>New</b>
<i>Perform Site</i>	CMC	CMC
<i>Method</i>	Turbidimetric	Turbidimetric
<i>PT Reference Range</i>	9.5-12.8 seconds	9.5-12.8 seconds
<i>Therapeutic INR Reference Range</i>	<b>0.83-1.13</b>	<b>0.85-1.14</b>
<i>Specimen</i>	Blue top- Sodium Citrate	Blue top- Sodium Citrate
<i>Volume Required</i>	4 mL	4 mL
<i>Specimen Stability</i>	24 hours	24 hours

  
 \_\_\_\_\_ Date 8/20/24  
 Laboratory Medical Director or Designee

  
 \_\_\_\_\_ Date 8/20/2024  
 Kristen Pellor, Core Laboratory Manager