

## Chemistry

Test Name	Stability	Container	Test Name	Stability	Container
Acetaminophen	72 R	SST	FSH	48 R	SST
Albumin	72 R	PST/SSST	FT3	48 R	SST
Alcohol	72 R	SST	FT4	48 R	SST
ALP	72 R	PST/SSST	Gentamicin	72 R	SST
ALT	72 R	PST/SSST	GGT	72 R	SST
Ammonia	3 ICE/30min T	Pink/Lav	Glucose	72 R	PST/SSST/Grey
Amylase	72 R	PST/SSST	HA1C	7 Days	Lavender
AST	72 R	PST/SSST	HDL	72 R	PST/SSST
Blood Gas	30 min T	Syringe	Hepatitis	48 R	SST
B12	24 R	SST	HIV	72 R	Dedicated Pink
BHCG	48 R	PST/SSST	Ion Calcium	30 min R	Green
BHCG reflex	48 R	Red	Insulin	24 R	SST
BNP	24 R	Pink/Lav	K	72 R	PST/SSST
BUN	72 R	PST/SSST	Lactic Acid	1 T	Grey
Calcium	72 R	PST/SSST	LDH	72 R	SST
Carbamazepine	72 R	SST	LDL	72 R	PST/SSST
Carbon Monoxide	7 days R T	Green/Syringe	LH	48 R	SST
CEA	48 R	SST	Lipase	72 R	PST/SSST
Cholesterol	72 R	PST/SSST	Lithium	72 R	SST
CK	12 R	PST/SSST	Lyme	7 Days	SST
CKMB	48 R	PST/SSST	Magnesium	72 R	PST/SSST
Cl	72 R	PST/SSST	Myoglobin	24 R	PST/SSST
CO2	72 R	PST/SSST	Na	72 R	PST/SSST
Cortisol	48 R	PST/SSST	Osmolality	72R / 24R	SST/Urine
Creatinine	72 R	PST/SSST	Phenobarbital	72 R	SST
CRP/CRPHS	72 R	PST/SSST	Phenytoin	72 R	SST
CSF Glucose	5 Days	CSF vial	Phosphorus	72 R	PST/SSST
CSF Protein	72 R	CSF vial	Prealbumin	72 R	SST
DBILI	72 R	PST/SSST	Procalcitonin	48 R	PST/SSST
Digoxin	72 R	SST	Progesterone	48 R	Red
Sens Estradiol	8 T / 72R	SST	Prolactin	48 R	SST
Fe (Iron)	72 R	SST	PSA	24 R	SST
Ferritin	48 R	SST	PTH intact	4T / 8R	SST
Folate	48 R	SST	PTH IO	8 T	Pink/Lav

## Hematology

Test Name	Stability	Container	Test Name	Stability	Container
BFCC	24 R	Specialty Tubes	QFTB	28 Days R	Specialty Tubes
BFCRY	24 R	SST	RF	72 R	SST
CBGD	24 R	SST	Rubella	48 R	SST
CSFCC	1 T	SST	Salicylate	72 R	SST
D Dimer	4 R	SST	Syphilis (RPR)	4 Days	SST
ESR	12R / 4T	PST/SSST	TBLI	72 R	PST/SSST
Malaria	24 R	SST	Testosterone	48 R	SST
Mono	48R SST/24R LAV	SST	TGAB	48 R	SST
Fibrinogen	4 R	SST	Theophylline	72 R	SST
PT/INR	24 R	SST	Total T3	48 R	SST
PTT	4 R	SST	Total T4	24 R	SST
Retic	24 R	PST/SSST	Total Protein	72 R	PST/SSST
Urinalysis	24R / 1T	SST	TPO	48 R	SST

## Blood Bank

Test Name	Stability	Container
DAT	24 R	SST
Fetal Screen	48 R	PST/SSST
Fetal Stain	24 R	SST
TS	72 R	SST

\*Stabilities in hours unless otherwise noted

Key	
R	Refrigerated
T	Room Temp

Updated 9/2023, 9/2024, 3/2025 KT

Refer to Cayuga Lab Test Catalog for additional container types and stability times.

MD 9/15/17

Reviewed by:  
Elizabeth Plocharczyk, MD

MAR 18 2025

X