

Preserved 24 Hour Urine Acceptability Chart

Test	Preservative/Acceptable (Yes/No)				
	Unpreserved	Acetic Acid	Boric Acid	Sodium Carbonate	6N HCL
Urine Creatinine LAB712	Yes	No*	Yes	No*	Yes
Urine & Serum Creatinine Clearance LAB1765**	Yes	No*	Yes	No*	Yes
Urine Total Protein LAB441	Yes	No	No*	No	No
Urine Microalbumin LAB410	Yes	No	No*	No	No
Urine Urea Nitrogen LAB449	Yes	Yes	Yes	Yes	Yes
Urine & Serum Urea Nitrogen Clearance LAB749**	Yes	Yes	Yes	Yes	Yes
Urine Electrolytes LAB19204	Yes	No*	No*	No	No
Urine Potassium LAB436	Yes	No*	No*	No*	No*
Urine Sodium LAB446	Yes	No*	No*	No	No*
Urine Chloride LAB375	Yes	No*	No*	No*	No

*If unable to perform in house testing but able to send out, cancel and order a miscellaneous test (GO23910). Refer to *Mayo Clinic Labs (MCL) Urine Preservative Chart* for test codes and additional preservative information.

**Tests require both serum and urine specimens