



Cayuga Medical Center Laboratory

Specimen Collection Supply Order Form

Please allow 3 to 5 business days for supplies to be filled.
 Please direct supply inquiries to Client Services (607) 274-4474
 Please fax supply requests to (607) 274-4481

Date: _____

Provider Office: _____ Contact Person: _____

Office Address: _____ Provider Phone Number: _____

| Stock # | ITEM | Quantity | Stock # | ITEM | Quantity |
|--|---|----------|--|--|----------|
| 001249 | 100/bag Biohazard Bags | | L00166 | 6pack Glucola 50gm | |
| | each Amber bags (protect from light) | | Urine Clinical Testing Specimen Collection | | |
| L01464 | 100/bag Vacutainer Needle Hubs | | 000754 | each "Nun's hat" toilet inserts | |
| 005076 | 50/box 23g Butterfly needle (lt blue) | | 000302 | 100/bag Orange Clinical sterile specimen cups | |
| 005077 | 50/box 21g Butterfly needle (green) | | L00359 | each 24 hour Urine jug | |
| L00119 | 48/box 22g Straight needle (black) | | 000305 | 50/box Blue cap vacutainer system urine collection kit | |
| L00121 | 48/box 21g Straight needle (green) | | L02214 | 100/flat Conical Urine Tubes (Yellow tube) | |
| 008203 | 25/box Tourniquet | | Microbiology Specimen Collection | | |
| Blood Collection Tubes: | | | 364953 | 50/box Urine C&S Tube/Transfer Straw kit | |
| 008425 | 100/ flat RED 6.0 mL | | White Culture swab (each) | | |
| 545070 | 100/ flat Lt BLUE (Na Citrate) 2.7 mL | | Red Culture swab (each) | | |
| 545071 | 100/ flat GOLD SST 5.0ml | | Green Culture swab (each) | | |
| 545074 | 100/ flat LAVENDER 4.0ml | | Orange Charcoal swab-Pertussis (each) | | |
| 545025 | 100/ flat PINK 6.0 mL | | Anaerobic (Port-A-Cult) swab (each) | | |
| 006758 | 100/flat CLEAR non-additive 6.0 ML tube | | M5 Viral Media swab (each) | | |
| 008424 | 100/ flat GRAY 6.0 mL | | BD Affirm VIII Transport kits (box) | | |
| Specialty Blood Collection Tubes | | | Aptima- GC/CHL URINE collection tubes (each) | | |
| each | Dark Blue/Purple EDTA (metal free) | | Aptima-Endocervical/Urethral swab (each) | | |
| each | Dark Blue/Red Serum tube (metal free) | | Stool Ova/Parasite kit (each) | | |
| each | Pearl Top EDTA plasma separator tube | | Cytology-Gynecological Specimen Collection | | |
| each | Green Sodium Heparin plasma tube | | 25/pack Wallach Pappette (cyto broom) | | |
| each | Quantiferon kit | | 25/pack Medscan CytoBrush | | |
| Pathology Specimen Collection | | | 25/flat ThinPrep PAP vials | | |
| Formalin in cup - 20ml | | | 1/each Soft ECC-Endocervical Curette | | |
| Formalin in cup - 30ml | | | 1/each Soft Biopsy-Gynecological Biopsy Device | | |
| Formalin in cup - 120 ml | | | Cytology- NON Gynecological Specimen Collection | | |
| Laboratory Test Requisitions | | | 20/pack CytoLyt Solution Vial | | |
| Laboratory GYN Cytology Requisition Form (13071) | | | Red Top Slide Jar (95% alcohol) | | |
| Laboratory Microbiology Requisition Form (13053) | | | Blue Frosted Microscope slide | | |
| Laboratory Surgical Pathology Form (13042) | | | Other Requests: | | |
| Laboratory Test General Requisition Form (13040) | | | | | |

Filled by: _____

Received Time/Date Stamp: _____