



A Member of Cayuga Health System

Laboratory Technical Bulletin

Date: 3/20/2025

Effective Date: 3/22/2025

Title: New Test Order: High Sensitivity Troponin-I, Random

Test: High Sensitivity Troponin-I

Explanation: A random High Sensitivity Troponin- I test will be orderable after the baseline and 1 hour troponin (Troponin I, Timed) have been ordered. The random test will not have a calculated delta or a critical flag.

	New
<i>Perform Site</i>	CMC Lab
<i>Test Code</i>	Troponin I (LAB747)
<i>Units</i>	pg/mL
<i>Reference Range</i>	Males <20 Females <15
<i>Reportable Range</i>	3-24,000
<i>Critical Values</i>	N/A
<i>Specimen</i>	Green PST only – invert 8-10x

Reviewed by:
Elizabeth Plocharczyk, MD
MAR 19 2025
x

Date _____
Laboratory Medical Director or Designee

Date 3/20/2025
Kristen Pellor, Core Laboratory Manager

Date 3/21/25
Michelle Walters, Laboratory Director