

General Requisition / Downtime Form

☐ STAT

☐ Routine

Name				Chart Label (preferred)	
DOB		Location/Rm. #/Unit			
Order Date & Time	Collected by	Collection Date & Time	Ordering Provider	Provider Signature	
Diagnosis codes				Results entered by	

For collection containers, refer to the test catalog for the testing facility.

<https://www.testmenu.com/arnothhealth> or <https://cayugamedlab.testcatalog.org/>

STAT TESTS (AVAILABLE DURING DOWNTIME) (C=Cayuga, A=Arnot)											
✓	Test Name	Specimen type	✓	Test Name	Specimen type	✓	Test Name	Specimen type	✓	Test Name	Specimen type
	Acetaminophen	SST		Carbon monoxide	Dark Green		Flu A/B	Swab		Salicylate	SST
	Activated partial thromboplastin time (PTT)	Lt. Blue		CBC, diff	Purple		Fungal cult/ smear	Swab		Stool occult blood	Card, cup, or Hemosure
	AFB cult/smear	Swab		CBC, no diff	Purple		HCG pregnancy	SST		Thyroid stimulating hormone (TSH)	SST
	Alcohol	SST		Comprehensive metabolic profile	SST/MINT		HIV 1&2	Pink		Troponin, baseline	Mint
	Ammonia	C-Pink A- Dark Green		COVID, Molecular	C- Aptima Orange		INR	Lt. Blue		Troponin, 1 hour	Mint
	Amylase	SST		COVID/Flu/RSV	Cepheid Vial		Ionized calcium	Dark Green		Troponin, 3 hour	Mint
	Arterial blood gas	Heparin syringe		CRP	SST		Lactic acid	Gray		Type & screen	
	B type natriuretic peptide (BNP)	Purple		CSF culture/smear	CSF Vial		Lipase	SST		If units are needed, use blood bank form	Pink
	Basic metabolic profile	SST		CSF glucose	CSF Vial		Magnesium	SST		Urinalysis with reflex culture	Cup, Urine collection tubes
	Blood culture	Blood culture bottles		CSF protein	CSF Vial		MRSA nasal PCR	Swab		Urine drug screen	Cup
	Body fluid, cell count	Syringe or sterile cup		D dimer quantitative	Lt. Blue		Rapid COVID	Swab		Venous blood gas	Heparin syringe or dark green
	Body fluid, crystals	Syringe or sterile cup		Direct bilirubin	SST		Rapid strep	Swab			
	Body fluid, glucose	Syringe or sterile cup		ESR/Sed rate	Purple		RSV	Swab			
	Body fluid, PH	Syringe or sterile cup		Fibrinogen	Lt. Blue		Rupture of membranes	Call Lab			
ROUTINE TESTS (DELAYED DURING DOWNTIME)											
	Estradiol, sensitive	C- SST A- SST/Mint		Hepatitis B surface antigen	SST		Progesterone	C- Red A- Red/SST		Thyroid peroxidase Ab	SST
	Ferritin	SST		Hepatitis C antibody	SST		Prolactin	C- SST A- SST/Mint		Transferrin	SST
	Fetal fibronectin	Call Lab		Insulin	SST		PSA Diagnostic	SST		Uric acid	SST
	Folate	C- SST A- SST/Mint		Iron	SST		PSA screening	SST		Valproic acid	SST
	Follicle Stimulating Hormone	C- SST A- SST/Mint		LDH	C- SST A- SST/Mint		PTH intact	C- SST A- SST/Mint		Vancomycin, peak	C- SST A- SST/Mint
	Free T3	C- SST A- SST/Mint		Lipid panel	C- SST A- SST/Mint		Reticulocyte count	Purple		Vancomycin, random	C- SST A- SST/Mint
	Free T4	C- SST A- SST/Mint		Lithium	SST					Vancomycin, trough	C- SST A- SST/Mint
	Gonorrhea/Chlamydia	Aptima Urine		Luteinizing hormone	C- SST A- SST/Mint		Salicylate	SST		Vit.D total 25(OH)	SST
	GGTP	C- SST A- SST/Mint		Osmolality, serum	SST		Syphilis IgG	SST		Vitamin B12	SST
	Hemoglobin A1C	Purple		Osmolality, urine	SST		T3, total	C- SST A- SST/Mint		Wound culture	Swab
	Hepatitis B surface antibody	SST		Prealbumin	C- SST A- SST/Mint		T4	C- SST A- SST/Mint		Source:	

Additional tests not listed(see test catalog for collection container):

Place downtime labels here

Document #

Version #

Date of Last Revision